

# Electroporation protocol for 3T3 cells

GTporator®

Transfection protocol

Protocol No. 09/2008-007

<b>Cell line</b>	3T3 (ATCC CCL-163)
<b>Washing solutions</b>	Phosphate buffered saline (PBS), pH 7.4, GTporator®-M
<b>Cell count</b>	1-3 x 10 <sup>6</sup>
<b>Electroporation solution</b>	GTporator®-M
<b>Cuvette</b>	2 mm gap width
<b>Volume</b>	80 µl
<b>Temperature</b>	Room temperature
<b>DNA</b>	5 µg in water

## Instrument settings

<u>Single amplitude</u>		<u>Radio-frequency</u>	
Protocol	Square-wave	Voltage	230 V
Voltage	200 V	Frequency	40 kHz
(Capacitance	960 µF)	Pulses	5 x 4 ms
Pulse time	20 ms	Interval	1 s
Pulse number	1	Modulation	100%

1. Trypsinize subconfluent cells
2. Wash the cells once in 1 ml PBS, and once in 300 µl of GTporator® solution
3. Resuspend the cells in 80 µl GTporator® solution with 5 µg of plasmid DNA
4. Transfer into a 2 mm electroporation cuvette
5. Electroporate
6. Immediately add pre-warmed culture medium with serum, transfer the cells into cell culture flask and place in incubator
7. Test the gene expression 24 h after electroporation

### Notes:

- The overall incubation time in the GTporator® solution must not exceed 20 minutes to guarantee a successful electroporation!
- The best results are obtained with DNA purified by anion-exchange chromatography.
- DNA should be dissolved in sterile water, using the TE buffer slightly lowers the survival rate.

### Typical results:

Survival rate in control cells 60%

Transfection rate: 60-75% with variations depending on the quality of DNA.

Results were measured 24 hours after transfection.

GTporator® is a registered trademark of the Protean s.r.o. company.

Protean s.r.o., Dolní 2, 370 04 České Budějovice, Czech Republic, Tel.: +420 724 731 831

**proteam**

E-mail: [info@protean.cz](mailto:info@protean.cz), web: [www.protean.cz](http://www.protean.cz)